

Product datasheet for **GE100057**

CRISPRa SAM Vector Kit, CRISPR/Cas9 Gene Activation

Product data:

Product Type: CRISPR Vectors

Function: CRISPRa SAM vector kit

Features:

- Robust CRISPR gene activation system
- Synergistic activation by three activation domains, VP64, p65 and HSF1
- Two vector system, dCas9-VP64-gRNA(MS2) (SKU GE100055) and MS2-p65-HSF1 (SKU GE100056)
- Scramble gRNA control (SKU GE100058) is included
- CF3 primer (GE100008) can be used to sequence target sequence after cloning
- Complement to gene overexpression with [cDNA clones](#)

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.



[View online »](#)