


Product datasheet for **GE100059**

pCas-Guide-CRISPRi gene Interference vector

Product data:

- Product Type:** CRISPR Vectors
- Function:** All-in-one CRISPRi vector
- Features:**
- All-in-one CRISPRi vector system
 - gRNA driven by U6 promoter (after target sequence cloned)
 - dCas9-KRAB-MeCP2 fusion protein expression by CMV promoter
 - CF3 primer (GE100008) can be used to sequence target sequence after cloning

Product images:



[View online >](#)