

Product datasheet for **GE100061**


pCas-Guide-Dead Mutant (D10A, H840A)

Product data:

Product Type:	CRISPR Vectors
Function:	Cas9 Nickase
Features:	All-in-one vector <ul style="list-style-type: none"> Expresses human codon-optimized Cas9 D10A/H840A mutants, C-terminal Myc-DDK tagged BamH I and BsmB I can be used for genomic target sequence cloning Expresses gRNA containing the inserted target sequence and the gRNA scaffold CF3 primer (GE100008) can be used to sequence the target sequence after cloning

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online >](#)