

Product datasheet for **GE100062**

pCas-Guide-Nickase (H840A)

Product data:

Product Type: CRISPR Vectors


Function: Cas9 Nickase

Features: All-in-one vector

- Expresses human codon-optimized Cas9 H840A mutants, C-terminal Myc-DDK tagged
- BamH I and BsmB I can be used for genomic target sequence cloning
- Expresses gRNA containing the inserted target sequence and the gRNA scaffold
- CF3 primer (GE100008) can be used to sequence the target sequence after cloning

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)