


Product datasheet for **GE100066**

pAAVS1-mRFP-BSD-DNR

Product data:

- Product Type:** CRISPR Vectors
- Function:** AAVS1 donor vector
- Features:** Donor vector for targeted transgene insertion into human AAVS1 site
- Contains left and right AAVS1 homologous sequences (550 bp each side)
 - Multiple cloning sites for transgene cloning, compatible for shuttling from the TrueORF clones
 - Blasticidin mammalian selection under PGK promoter
 - Targeted insertion when cotransfected with pCas-Guide-AAVS1 (GE100023)

Product images:



[View online >](#)

Kozac
Consensus

EcoR I Sal I Kpn I RBS Sgf I

CTATAGGGCGGCCGGGAATTCTGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGGATCGCC

Asc I Rsr II Mlu I Not I

GGCGGCCAGATCTCAAGCTTAAGTACTAGCTAGCGGACCGACGGGTACGGGCCGCTC ATG GTG TCT

M V S

mRFP tag Pme I Fse I

--- --- GGA CTC AGA GTT TAA ACGGCCGGCC

G L R V stop