


Product datasheet for **GE100067**

pRosa26-tGFP-Puro-DNR

Product data:

Product Type:	CRISPR Vectors
Function:	ROSA26 donor vector
Features:	<ul style="list-style-type: none"> • Donor vector for targeting tGFP into Mouse ROSA26 locus • Contains left and right ROSA26 homologous sequences (550 bp each side) • Puromycin mammalian selection under PGK promoter • Targeted insertion when cotransfect with pCas-Guide-ROSA26 (GE100050) (SKU GE100050)
Disclaimer:	These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)

EcoRI
BamHI KpnI
RBS
Kozak
Consensus
SgfI

CTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGGATCGCC

AscI
MluI
NotI
GGCGGCCAGATCTCAAGCTTAAGTGTAGCGGACCGACGCGT ACG CGG CCG CTC GAG ATG
T R P L E M

tGFP Tag
PmeI
GAG AGC GAC --- --- GGT GAA GAA AGA GTT TAA AC

E S D - - - G E E R V Stop