



## Product datasheet for **GE100070**

### pAAVS1-EF1a-tRFP-Puro-DNR

#### Product data:

- Product Type:** CRISPR Vectors
- Function:** AAVS1 donor vector
- Features:**
- Donor vector for targeted transgene insertion into human AAVS1 site
  - Contains left and right AAVS1 homologous sequences (550 bp each side)
  - Multiple cloning sites for transgene cloning, compatible for shuttling from the TrueORF clones
  - Puromycin mammalian selection under PGK promoter
  - Targeted insertion when cotransfected with pCas-Guide-AAVS1 (SKU GE100023)

#### Product images:



View online »