

Product datasheet for **GE100075**


pLenti-Cas9-Myc-DDK-P2A-Puro

Product data:

- Product Type:** CRISPR Vectors
- Function:** Cas9 only vectors
- Features:** Cas9 only vector with puro selection in lenti backbone
- 3rd generation Lenti vector, ideal for hard-to-transfect cells
 - Cas9-Myc-DDK fusion protein under CMV promoter for mammalian expression after transfection or transduction
 - P2A-Puro is within the Lentiviral packaging cassette and would allow Puromycin selection in both the transfected and the transduced cells.
 - Lenti particles can be produced using [high efficiency lenti packaging kits](#)

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



[View online »](#)