

Product datasheet for **GE100085**

pCas-Guide-GFP-CRISPRi gene interference vector


Product data:

Product Type: CRISPR Vectors

- Features:
- All-in-one CRISPRi vector system
 - gRNA driven by U6 promoter (after target sequence cloned)
 - dCas9-KRAB-MeCP2 fusion protein expression by CMV promoter
 - CF3 primer (GE100008) can be used to sequence target sequence after cloning
 - CF3 primer (GE100008) can be used to sequence target sequence after cloning
 - Contains tGFP reporter gene under EF1a promoter

Disclaimer: These products are manufactured and supplied by OriGene under license from ERS.

Product images:



Plasmid map of pCas-Guide-GFP-CRISPRi gene interference vector (GE100085)



[View online >](#)