

Product datasheet for **GE100088**

Sequencing Primer, EF51-F

Product data:

Product Type:	DNA Sequencing Primers
Description:	Forward primer to sequence targets cloned in expression vectors.
Forward Sequence:	5' CAAGCCTCAGACAGTGGTTC 3'
Applicable Vectors:	PS100072, PS100084, PS100085, PS100106, PS100120, PS100121, PS100132
Components:	1 vial of lyophilized sequencing primer (1 nmol, sufficient for 100 sequencing reactions)
Quality Control:	The primer has been tested to generate satisfactory sequencing data on ABI 3730.
Storage:	Store at -20°C.
Stability:	The lyophilized and suspended primer is stable for one year from date of shipping when stored at -20°C



[View online »](#)