

Product datasheet for **GE100090**

Sequencing Primer, **GE100002-F**

Product data:

| | |
|---------------------|--|
| Product Type: | DNA Sequencing Primers |
| Description: | Forward primer to sequence targets cloned in expression vectors. |
| Forward Sequence: | 5' GATACAAGGCTGTTAGAGAG 3' |
| Applicable Vectors: | PS100002 |
| Components: | 1 vial of lyophilized sequencing primer (1 nmol, sufficient for 100 sequencing reactions) |
| Quality Control: | The primer has been tested to generate satisfactory sequencing data on ABI 3730. |
| Storage: | Store at -20°C. |
| Stability: | The lyophilized and suspended primer is stable for one year from date of shipping when stored at -20°C |



[View online »](#)