

Product datasheet for **GE100091**

Sequencing Primer, **GE100010-F**

Product data:

Product Type: DNA Sequencing Primers

Description: Forward primer to sequence targets cloned in expression vectors.

Forward Sequence: 5' CATGATTCCTTCATATTTGC 3'

Components: 1 vial of lyophilized sequencing primer (1 nmol, sufficient for 100 sequencing reactions)

Quality Control: The primer has been tested to generate satisfactory sequencing data on ABI 3730.

Storage: Store at -20°C.

Stability: The lyophilized and suspended primer is stable for one year from date of shipping when stored at -20°C



[View online »](#)