

## Product datasheet for **GE100104**

### Sequencing Primer, PR100001-F

#### Product data:

Product Type:	DNA Sequencing Primers
Description:	Forward primer to sequence targets cloned in expression vectors.
Forward Sequence:	5' GGCTATTGGCCATTGCATAC 3'
Applicable Vectors:	PR100001, PR100002
Components:	1 vial of lyophilized sequencing primer (1 nmol, sufficient for 100 sequencing reactions)
Quality Control:	The primer has been tested to generate satisfactory sequencing data on ABI 3730.
Storage:	Store at -20°C.
Stability:	The lyophilized and suspended primer is stable for one year from date of shipping when stored at -20°C



[View online »](#)