

Product datasheet for GE100119

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Sequencing Primer, GST-R

Product data:

Product Type: DNA Sequencing Primers

Description: Reverse primer to sequence targets cloned in expression vectors.

5' CCATTTATCACCTTCATCGC 3' **Forward Sequence:**

Applicable Vectors: PS100126

1 vial of lyophilized squencing primer (1 nmol, sufficient for 100 sequencing reactions) Components:

Quality Control: The primer has been tested to generate satisfactory squencing data on ABI 3730.

Store at -20°C. Storage:

Stability: The lyophilized and suspended primer is stable for one year from date of shipping when

stored at -20°C