

Product datasheet for **GE100125**

Sequencing Primer, mGFP-R

Product data:

Product Type:	DNA Sequencing Primers
Description:	Reverse primer to sequence targets cloned in expression vectors.
Forward Sequence:	5' GAAGCACTGGATGCCGTAGC 3'
Applicable Vectors:	PS100040, PS100071, PS100084, PS100093, PS100102, PS100113, PS100119, PS100121
Quality Control:	The primer has been tested to generate satisfactory sequencing data on ABI 3730.
Storage:	Store at -20°C.
Stability:	The lyophilized and suspended primer is stable for one year from date of shipping when stored at -20°C



[View online »](#)