

## Product datasheet for **GE100128**

### Sequencing Primer, P2A-DDK-R

#### Product data:

Product Type:	DNA Sequencing Primers
Description:	Reverse primer to sequence targets cloned in expression vectors.
Forward Sequence:	5' GGACGTCGTCGCGGGTGGCG 3'
Applicable Vectors:	PS100092, PS100101, PS100109, PS100110, PS100112, PS100114, PS100118, PS100120
Components:	1 vial of lyophilized sequencing primer (1 nmol, sufficient for 100 sequencing reactions)
Quality Control:	The primer has been tested to generate satisfactory sequencing data on ABI 3730.
Storage:	Store at -20°C.
Stability:	The lyophilized and suspended primer is stable for one year from date of shipping when stored at -20°C



[View online »](#)