

## Product datasheet for **GE100137**

### Sequencing Primer, XL39-R

#### Product data:

**Product Type:** DNA Sequencing Primers

**Description:** Reverse primer to sequence targets cloned in expression vectors.

**Forward Sequence:** 5' ATTAGGACAAGGCTGGTGGG 3'

**Applicable Vectors:** PS100001, PS100002, PS100004, PS100005, PS100006, PS100007, PS100008, PS100009, PS100011, PS100012, PS100013, PS100014, PS100016, PS100017, PS100018, PS100019, PS100020, PS100022, PS100024, PS100025, PS100026, PS100033, PS100035, PS100047, PS100048, PS100049, PS100051, PS100052, PS100055, PS100056, PS100057, PS100058, PS100061, PS100066, PS100067, PS100089, PS100111

**Components:** 1 vial of lyophilized sequencing primer (1 nmol, sufficient for 100 sequencing reactions)

**Quality Control:** The primer has been tested to generate satisfactory sequencing data on ABI 3730.

**Storage:** Store at -20°C.

**Stability:** The lyophilized and suspended primer is stable for one year from date of shipping when stored at -20°C



[View online »](#)