

Product datasheet for **GP10108-100**

P2rx3 Guinea Pig Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Rat
Host:	Guinea Pig
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	purinergic receptor P2X 3
Database Link:	Entrez Gene 81739 Rat
Synonyms:	P2X purinoceptor 3, ATP receptor, Purinergic receptor



[View online »](#)