

Product datasheet for **GP14100-50**

Trpv1 (C-term) Guinea Pig Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Rat
Host:	Guinea Pig
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	transient receptor potential cation channel, subfamily V, member 1
Database Link:	Entrez Gene 83810 Rat
Synonyms:	VR1, Capsaicin receptor, OTRPC1



[View online »](#)