

## Product datasheet for **GP14103-150**

### **TAC1 (1-11) Guinea Pig Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Human, Mouse, Rat
Host:	Guinea Pig
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens tachykinin precursor 1 (TAC1), transcript variant beta
Database Link:	<a href="#">Entrez Gene 21333 Mouse</a> <a href="#">Entrez Gene 24806 Rat</a> <a href="#">Entrez Gene 6863 Human</a>
Synonyms:	Substance-P, SP
Protein Families:	Druggable Genome, Secreted Protein

#### **Product images:**



