

Product datasheet for **GP14103-50**

TAC1 (1-11) Guinea Pig Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Reactivity:	Human, Mouse, Rat
Host:	Guinea Pig
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	Homo sapiens tachykinin precursor 1 (TAC1), transcript variant beta
Database Link:	Entrez Gene 21333 Mouse Entrez Gene 24806 Rat Entrez Gene 6863 Human
Synonyms:	Substance-P
Protein Families:	Druggable Genome, Secreted Protein

Product images:



