


## Product datasheet for **GP14104-50**

### **PGP9.5 (UCHL1) (175-191) Guinea Pig Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Human, Mouse, Rat
Host:	Guinea Pig
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	ubiquitin C-terminal hydrolase L1
Database Link:	<a href="#">Entrez Gene 7345 Human</a>
Synonyms:	UCH-L1, PGP 9.5, Ubiquitin thioesterase L1, Neuron cytoplasmic protein 9.5

#### **Product images:**



[View online »](#)

