


Product datasheet for **GP14108-100**

TRPM6 Guinea Pig Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Host:	Guinea Pig
Clonality:	Polyclonal
Formulation:	State: Serum
Gene Name:	transient receptor potential cation channel subfamily M member 6
Synonyms:	Channel kinase 2


Product images:



TRPM6 staining cryostat section of rat kidney. Tissues were stained with using Cy3-conjugated secondary antibodies.



[View online »](#)



TRPM6 staining cryostat section of rat kidney. Tissues were stained with using Cy3-conjugated secondary antibodies.