


Product datasheet for **GP14110-100**

TAC1 (1-11) Guinea Pig Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC
Host:	Guinea Pig
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens tachykinin precursor 1 (TAC1), transcript variant beta
Database Link:	Entrez Gene 6863 Human
Synonyms:	Substance-P
Protein Families:	Druggable Genome, Secreted Protein


Product images:




Immuofluorescent detection of Affinity Purified Substance P in rat spinal cord dorsal horn (red fluorescence).




[View online »](#)




Immuofluorescent detection of Affinity Purified Substance P in rat DRG (red fluorescence).




SP staining of rat spinal cord dorsal horn.
Neuroscience Letters Volume 501, Issue 1, 21
August 2011, Pages 4-9.



Immuofluorescent detection of Affinity Purified Substance P in rat spinal cord dorsal horn (red fluorescence).



Immuofluorescent detection of Affinity Purified Substance P in rat DRG (red fluorescence).



SP staining of rat spinal cord dorsal horn.
Neuroscience Letters Volume 501, Issue 1, 21
August 2011, Pages 4-9.