


Product datasheet for **GT41001-100**

ALS2 (C-term) Goat Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Reactivity:	Human
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens ALS2, alsin Rho guanine nucleotide exchange factor (ALS2), transcript variant 2
Database Link:	Entrez Gene 57679 Human
Synonyms:	ALS2CR6, KIAA1563
Protein Families:	Druggable Genome
Protein Pathways:	Amyotrophic lateral sclerosis (ALS)


Product images:



Alsin (dilution: 3.8ug/ml)



[View online »](#)



Alsin (dilution: 3.8ug/ml)