


## Product datasheet for **GT41011-100**

### **GABA A Receptor beta 3 (GABRB3) Goat Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Reactivity:	Human
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens gamma-aminobutyric acid type A receptor beta3 subunit (GABRB3), transcript variant 1
Database Link:	<a href="#">Entrez Gene 2562 Human</a>
Synonyms:	GABRB-3, GABA A receptor subunit beta-3
Protein Families:	Druggable Genome, Ion Channels: Cys-loop Receptors, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction


#### **Product images:**



GABA A Receptor beta 3 (dilution: 2.5ug/ml) staining of paraffin embedded Human Cortex. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.



[View online »](#)



GABA A Receptor beta 3 (dilution: 2.5ug/ml) staining of paraffin embedded Human Cortex. Steamed antigen retrieval with citrate buffer pH 6, AP-staining.