


## Product datasheet for **GT41015-100**

### **MSX1 (N-term) Goat Polyclonal Antibody**

#### **Product data:**

Product Type:	Primary Antibodies
Applications:	ELISA, IHC
Reactivity:	Human
Host:	Goat
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	Homo sapiens msh homeobox 1 (MSX1)
Database Link:	<a href="#">Entrez Gene 4487 Human</a>
Synonyms:	MSX-1, HOX7, Hox-7
Protein Families:	Druggable Genome, Transcription Factors


#### **Product images:**



MSX1 staining of staining of paraffin embedded Human Prostate (dilution: 2.5ug/ml). Steamed antigen retrieval with citrate buffer pH 6, AP-staining. Protocol on data-sheet.



[View online »](#)



MSX1 staining of staining of paraffin embedded Human Prostate (dilution: 2.5ug/ml). Steamed antigen retrieval with citrate buffer pH 6, AP-staining. Protocol on data-sheet.