

Product datasheet for **HCRT103**

TissueScan, Colon Cancer cDNA Array III

Product data:

Product Type:	TissueScan
Array Category:	Cancer Specific
Clinical Info:	/assets/documents/tissuescan/hcrt103.xls
Diagnosis Category:	Colon Cancer
Tissue:	Colon
Array Name:	Colon Cancer cDNA Array III


Product images:



TissueScan product image.




[View online »](#)



A Human Matched Colon Cancer real time PCR Gene Expression panel used in real-time PCR. The beta-actin contents in the 48 human cDNAs were detected using SYBR Green I protocol and performed in PE7700.

Fluorescence reading - cycle number.



A Human Matched Colon Cancer real time PCR Gene Expression panel was used in SYBR Green I based real-time PCR. The expression profile of the human Survivin transcript was detected using SYBR Green I protocol and performed in PE7700.

Normalized relative expression among 48 matched colon samples