

## Product datasheet for **HCRT501**

### TissueScan, Colon Cancer cDNA Array I

#### Product data:

Product Type:	TissueScan
Array Category:	Cancer Specific
Clinical Info:	/assets/documents/tissuescan/hcrt501.xls
Diagnosis:	Colon Cancer
Tissue:	Colon
Array Name:	Colon Cancer cDNA Array I


#### Product images:



TissueScan product image.




[View online »](#)



A Human Colon Cancer real time PCR Gene Expression panel used in real-time PCR. The beta-actin contents in the 48 human cDNAs were detected using SYBR Green I protocol and performed in PE7700.

Fluorescence reading - cycle number.



A Human Colon Cancer real time PCR Gene Expression panel was used in SYBR Green I based real-time PCR. The expression profile of the human Carbonic Anhydrase II (ACC#: NM\_000067) was detected using SYBR Green I protocol and performed in PE7700.

Normalized relative expression among 48 colon samples