

## Product datasheet for **HK200013**

### AAK1 Human qPCR Template Standard (BC002695)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC002695
Synonyms:	adaptor-associated kinase 1; AP2 associated kinase 1; DKFZp686F03202; DKFZp686K16132; FLJ23712; FLJ25931; FLJ31060; FLJ42882; FLJ45252; KIAA1048; MGC138170; MGC164568; MGC164570
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP211249) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)