

Product datasheet for **HK200015**

Acyl CoA Thioesterase 9 (ACOT9) Human qPCR Template Standard (BC012573)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC012573
Synonyms:	ACATE2; acyl-CoA thioesterase 9; acyl-Coenzyme A thioesterase 2, mitochondrial; CGI-16; mitochondrial Acyl-CoA Thioesterase; MT-ACT48; OTTHUMP00000023065
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP233391) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)