

Product datasheet for **HK200017**

ADAT2 Human qPCR Template Standard (BC037955)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC037955
Synonyms:	adenosine deaminase, tRNA-specific 2, TAD2 homolog (<i>S. cerevisiae</i>); DEADC1; deaminase domain containing 1; dj20N2; dj20N2.1; DKFZp686L1118; FLJ44213; OTTHUMP00000017329; OTTHUMP00000017330; TAD2; tRNA-specific adenosine deaminase 2 homolog
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP219060) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)