

Product datasheet for **HK200018**

Adenosine A3 Receptor (ADORA3) Human qPCR Template Standard (BC029831)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC029831
Synonyms:	A3AR; AD026; adenosine A3 receptor; bA552M11.5; OTTHUMP00000013476; OTTHUMP00000013739; OTTHUMP00000082917; RP11-552M11.7
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP231209) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)