

## Product datasheet for **HK200021**

### ALG10 Human qPCR Template Standard (BC070347)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC070347
Synonyms:	alpha2-glucosyltransferase; asparagine-linked glycosylation 10 homolog (yeast, alpha-1,2-glucosyltransferase); asparagine-linked glycosylation 10, alpha-1,2-glucosyltransferase homolog (S. p; derepression of ITR1 expression 2 homolog; DIE2; FLJ14751; KCR1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP216196) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)