

## Product datasheet for **HK200022**

### **ANAPC13 Human qPCR Template Standard (BC005398)**

#### **Product data:**

|                         |   |
|-------------------------|---|
| <b>Product Type:</b>    | qPCR Template Standards   |
| <b>ACCN:</b>            | BC005398  |
| <b>Synonyms:</b>        | anaphase promoting complex subunit 13; APC13; cyclosome subunit 13; DKFZp566D193; homolog of yeast Swm1; SWM1   |
| <b>Components:</b>      | 1 vial of lyophilized SybGREEN qPCR primer (HP211603) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| <b>Quality Control:</b> | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| <b>Storage:</b>         | Store at -20°C.   |
| <b>Stability:</b>       | The primer mix and template standard are stable for one year from date of shipping.   |



[View online »](#)