

Product datasheet for **HK200029**

Apg12 (ATG12) Human qPCR Template Standard (BC012266)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC012266
Synonyms:	APG12; Apg12 (autophagy, yeast) homolog; APG12 autophagy 12-like; APG12L; ATG12 autophagy related 12 homolog (<i>S. cerevisiae</i>); FBR93; HAPG12
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207980) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)