

## Product datasheet for **HK200031**

### ATP5L Human qPCR Template Standard (BC015128)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC015128
Synonyms:	ATP5JG; ATPase subunit G; ATP synthase, H <sup>+</sup> transporting, mitochond; ATP synthase, H <sup>+</sup> transporting, mitochondrial F1F0, subunit g; ATP synthase g chain, mitochondrial; F1F0-type ATP synthase subunit g; F1Fo-ATP synthase complex Fo membrane domain g subunit
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP209419) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)