

## Product datasheet for **HK200052**

### **C17orf81 (ELP5) Human qPCR Template Standard (AB013910)**

#### **Product data:**

Product Type:	qPCR Template Standards
ACCN:	AB013910
Synonyms:	C17orf81; DERP6; HSPC002; MST071; MSTP071
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP211582) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)