

## Product datasheet for **HK200083**

### **CARD14 Human qPCR Template Standard (BC018142)**

#### **Product data:**

<b>Product Type:</b>	qPCR Template Standards
<b>ACCN:</b>	BC018142
<b>Synonyms:</b>	bcl10-interacting maguk protein 2; BIMP2; CARD-containing MAGUK 2 protein; card-maguk protein 2; CARMA2; caspase recruitment domain family, member 14; caspase recruitment domain protein 14
<b>Components:</b>	1 vial of lyophilized SybGREEN qPCR primer (HP214753) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
<b>Quality Control:</b>	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
<b>Storage:</b>	Store at -20°C.
<b>Stability:</b>	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)