

## Product datasheet for **HK200085**

### CBR3 Human qPCR Template Standard (BC002812)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC002812
Synonyms:	carbonyl reductase (NADPH) 3; carbonyl reductase 3; hCBR3; NADPH-dependent carbonyl reductase 3; SDR21C2; short chain dehydrogenase/reductase family 21C, member 2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205162) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)