

Product datasheet for **HK200093**

MCP4 (CCL13) Human qPCR Template Standard (BC008621)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC008621
Synonyms:	chemokine (C-C; CK-beta-10; CKb10; MCP-4; MGC17134; monocyte chemoattractant protein 4; monocyte chemotactic protein 4; NCC-1; NCC1; new CC chemokine 1; SCYA13; SCYL1; small inducible cytokine A13; small inducible cytokine subfamily A (Cys-Cys), member 13
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP208545) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)