

## Product datasheet for **HK200094**

### LD78 beta (CCL3L1) Human qPCR Template Standard (BC027888)

#### Product data:

Product Type:	qPCR Template Standards
ACCN:	BC027888
Synonyms:	chemokine (C-C motif) ligand 3-like 1; D17S1718; G0S19-2; LD78; LD78BETA; MGC12815; MGC104178; MGC182017; MIP1AP; SCYA3L; SCYA3L1; small inducible cytokine A3-like 1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP213988) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 <sup>7</sup> copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)