

Product datasheet for **HK200095**

CD1 (CD1A) Human qPCR Template Standard (BC031645)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC031645
Synonyms:	CD1; CD1A antigen; CD1A antigen, a polypeptide; cluster of differentiation 1 A; cortical thymocyte antigen CD1A; differentiation antigen CD1-alpha-3; FCB6; HTA1; OTTHUMP00000018907; R4; T-cell surface antigen T6/Leu-6; T-cell surface glycoprotein CD1a; T6
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205560) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)