

Product datasheet for **HK200096**

CD48 Human qPCR Template Standard (BC030224)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC030224
Synonyms:	BCM1; BLAST; BLAST1; CD48 antigen (B-cell membrane protein); CD48 molecule; hCD48; mCD48; MEM-102; OTTHUMP00000025680; OTTHUMP00000060268; SLAMF2
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP205570) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)