

Product datasheet for **HK200097**

p16INK4A (CDKN2A) Human qPCR Template Standard (BC015960)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC015960
Synonyms:	ARF; CDK4 inhibitor p16-INK4; CDK4I; CDKN2; cell cycle negative regulator beta; CMM2; cyclin-dependent ki; cyclin-dependent kinase inhibitor 2A; INK4; INK4a; MLM; MTS1; OTTHUMP00000021147; OTTHUMP00000021148; p14; p14ARF; p16; p16INK4; p16INK4a; p19; TP16
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP226191) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)