

Product datasheet for **HK200099**

Centrin 1 (CETN1) Human qPCR Template Standard (BC029515)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC029515
Synonyms:	calcium binding protein; CEN1; centrin, EF-hand protein, 1; Centrin-1; centrin 1; CETN; EF-hand protein
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207437) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)