

Product datasheet for **HK200113**

CLACP (COL25A1) Human qPCR Template Standard (BC036669)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC036669
Synonyms:	CLAC; CLACP; collagen, type XXV, alpha 1; collagen-like Alzheimer amyloid plaque component; collagenous Alzheimer amyloid plaque component; OTTHUMP00000196081; OTTHUMP00000196082
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP230701) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)