

## Product datasheet for **HK200116**

### **INT11 (CPSF3L) Human qPCR Template Standard (AL136813)**

#### **Product data:**

|                         |   |
|-------------------------|---|
| <b>Product Type:</b>    | qPCR Template Standards   |
| <b>ACCN:</b>            | AL136813  |
| <b>Synonyms:</b>        | CPSF73L; INT11; INTS11; RC-68; RC68   |
| <b>Components:</b>      | 1 vial of lyophilized SybGREEN qPCR primer (HP212660) mix (1 nmol each primer, sufficient for 200 reactions)<br>1 vial lyophilized Gene Specific qPCR template standard: $50 \times 10^7$ copies, double-stranded DNA<br>1 vial containing 1 mL 10x PCR template standard dilution buffer |
| <b>Quality Control:</b> | The kit has been tested to generate satisfactory qPCR data on ABI 7900HT  |
| <b>Storage:</b>         | Store at -20°C.   |
| <b>Stability:</b>       | The primer mix and template standard are stable for one year from date of shipping.   |



[View online »](#)