

Product datasheet for **HK200125**

DBF4B Human qPCR Template Standard (BC016158)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC016158
Synonyms:	activator of S-phase kinase-like protein 1; ASKL1; chifb; chiffon homolog B; DBF4 homolog B; DBF4 homolog B (<i>S. cerevisiae</i>); Dbf4-related factor 1; DRF1; FLJ13087; MGC15009; ZDBF1B; zinc finger, DBF-type containing 1B
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP233133) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)