

Product datasheet for **HK200138**

EBAG9 Human qPCR Template Standard (BC017729)

Product data:

Product Type:	qPCR Template Standards
ACCN:	BC017729
Synonyms:	cancer associated surface antigen; EB9; estrogen receptor binding site associated, antigen, 9; estrogen receptor binding site associated antigen 9; PDAF; RCAS1
Components:	1 vial of lyophilized SybGREEN qPCR primer (HP207565) mix (1 nmol each primer, sufficient for 200 reactions) 1 vial lyophilized Gene Specific qPCR template standard: 50 X 10 ⁷ copies, double-stranded DNA 1 vial containing 1 mL 10x PCR template standard dilution buffer
Quality Control:	The kit has been tested to generate satisfactory qPCR data on ABI 7900HT
Storage:	Store at -20°C.
Stability:	The primer mix and template standard are stable for one year from date of shipping.



[View online »](#)